

**E006: Temperature monitoring devices**

<b>PQS code:</b>	<b>E006/092</b>
<b>Description:</b>	Irreversible Freeze Indicator
<b>Manufacturer's ref:</b>	LM Pro IN031
<b>Manufactured in:</b>	India
<b>Company:</b>	G-Tek Corporation Private Limited
<b>Address:</b>	3, MAHAVIR INDUSTRIAL ESTATE NEAR JALARAM MANDIR, Karelibagh, Vadodara - 390018 Gujarat, India
<b>Telephone:</b>	+91 98 242 61702
<b>Email:</b>	<a href="mailto:asu@gtek-india.com">asu@gtek-india.com</a>
<b>Web address:</b>	<a href="http://www.gtek-india.com">www.gtek-india.com</a>

**Specifications**

<b>Product type:</b>	Irreversible freeze indicator	<b>Display medium:</b>	LCD
<b>Mode of operation:</b>	Electronic	<b>Weight</b>	0,030 kg
<b>Minimum temperature:</b>	- 20°C	<b>Main material:</b>	Polycarbonate plastic
<b>Maximum temperature:</b>	+50°C	<b>Dimensions:</b>	7.0 L x 3.5 W x 1.3 H cm

**MECHANICAL AND ELECTRONIC DEVICES ONLY**

<b>Accuracy:</b>	± 0.5 °C	<b>Min logging interval:</b>	
<b>Power source:</b>	Non-replaceable battery	<b>Alarm type:</b>	Not applicable
<b>Battery type:</b>	Lithium Coin Cell Battery	<b>User interface:</b>	LCD
<b>Battery shelf life:</b>	3.5 years	<b>Programmability:</b>	Not applicable
<b>Activated life:</b>	42 months	<b>IP rating:</b>	
<b>Memory capacity:</b>			

**Comments:**

- 1) Device is delivered in sleep mode.
- (2) Activate the device by long press of key for 10 seconds.
- (3) Ensure the display Shows a check mark/tick mark before placing it with vaccine
- (4) If the device is exposed to temperature of -0.5 °C, ± 0.5 °C for 60 minutes ± 5 minutes (or as specified for the TYPE A/B/C), alarm is activated and display shows "X". This is irreversible.
- (5) For more details refer to the Quick start guide.

<b>Minimum order:</b>	100	<b>Year base price:</b>	2023
<b>Quality standard:</b>	- ISO 9001:2015 -	<b>Pieces per carton:</b>	
<b>Price per unit:</b>	\$6 for 100 units	<b>Volume per carton:</b>	0.0072 m3
		<b>Weight per carton:</b>	5 kg
		<b>Incoterms:</b>	EXW
<b>Verification report:</b>		<b>Verification laboratory:</b>	
<b>Current PQS status:</b>	<i>pre-qualified ; 15 May 2023</i>	<b>Valid until:</b>	<i>May 2026</i>

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.